



MICAP

## MICAP-OES 1000

Direct analysis of trace elements in Middle Distillate Fuels with a compact N<sub>2</sub> ICP-OES system

### Introduction

Middle distillates are a category of fuels produced from the distillation of crude oil. This includes products commonly referred to as diesel, marine diesel oil, kerosene and jet fuel. These colorless to light yellow fuels each differ slightly in their components based on the application. Elemental impurities present can impact critical properties such as combustion, corrosion, shelf life, and production of engine deposits. Inductively Coupled Plasma Optical Emission Spectroscopy (ICP-OES) is commonly utilized to monitor the trace elemental content of these various fuels.

This note describes the use of the simultaneous MICAP™ OES-1000 N<sub>2</sub> ICP-OES system for routine analysis of middle distillate fuels following the procedures outlined in ASTM Method D7111. Analytical performance of the MICAP is presented to demonstrate accurate and stable performance during the direct analysis of several fuel samples.

**RADOM**  
INSTRUMENTS

## Instrumentation

The ASTM standard test method D7111<sup>1</sup> covers the middle distillate products, which includes distillation fractions between 150 and 390 °C. The samples are analyzed directly against organometallic standards prepared in a kerosene solvent matrix. This requires the Radom Instruments MICAP-OES 1000 to be configured with an organic solvent introduction system. A double-pass cyclonic spray chamber and high-solids V-Groove nebulizer were connected to the standard one-piece torch for introduction of samples into the N<sub>2</sub> plasma.

Microwave energy in the MICAP is coupled into the Cerawave™ ring in a highly efficient process that creates the magnetic fields required to inductively couple energy into the robust N<sub>2</sub> plasma source. The 4 MP sCMOS camera simultaneously collects the resulting emission lines from the high-resolution spectrometer.

## Experimental Conditions

The analyses were performed on the simultaneous MICAP (Figure 1) utilizing the conditions listed in Table 1. Samples were introduced with an autosampler, with each sample analysis taking approximately 3.5 minutes to complete.

The MICAP in this work was uniquely configured with air (rather than the standard N<sub>2</sub>) as the nebulizer gas. This modification reduces the continuum background created by the combustion of carbon, resulting in significantly improved detection limits for several elements. Table 2 provides specifics on the analyte and Yttrium (Y) internal standard (I/S) wavelengths used for this analysis, as well as the detection limits obtained from 3x the standard deviation of 10 blank measurements.



Figure 1. MICAP-OES 1000 N<sub>2</sub> ICP-OES

Table 1. Operational conditions

| Parameter          | Value                          |
|--------------------|--------------------------------|
| Torch              | Quartz 1-piece, 1.5mm injector |
| Spray Chamber      | Double pass cyclonic           |
| Nebulizer          | V-Groove glass                 |
| Sample Tubing      | Solvaflex Blk/Blk (0.76 mm ID) |
| Drain Tubing       | Solvaflex Blu/Yel (1.52 mm ID) |
| Rinse              | PremiSolv™                     |
| Stabilization Time | 15 sec                         |
| Power              | 1000 W                         |
| Coolant Gas Flow   | 16 L/min N <sub>2</sub>        |
| Auxiliary Gas Flow | 0.4 L/min N <sub>2</sub>       |
| Nebulizer Gas Flow | 0.55 L/min Air                 |
| Peristaltic Pump   | 90 rpm                         |
| Plasma Viewing     | Axial                          |
| Camera Exposure    | 30 sec (125 ms @ 240 exps)     |
| # of Repeats       | 3                              |

Table 2. Analyte wavelengths and fuel detection limits

| Analyte Wavelength (nm) | Detection Limit (mg/Kg, ppm) |
|-------------------------|------------------------------|
| Ag 328.068              | 0.001                        |
| Al 396.152              | 0.007                        |
| B 249.677               | 0.022                        |
| Ba 585.368              | 0.018                        |
| Ca 393.366              | 0.0007                       |
| Cd 228.802              | 0.007                        |
| Cr 428.973              | 0.003                        |
| Cu 324.754              | 0.0007                       |
| Fe 259.940              | 0.006                        |
| K 766.490               | 0.020                        |
| Mg 279.553              | 0.0003                       |
| Mn 257.610              | 0.002                        |
| Mo 317.034              | 0.007                        |
| Na 588.995              | 0.002                        |
| Ni 300.249              | 0.016                        |
| P 213.618               | 0.380                        |
| Pb 283.305              | 0.035                        |
| Si 251.611              | 0.003                        |
| Sn 283.998              | 0.019                        |
| Ti 334.940              | 0.002                        |
| V 309.310               | 0.007                        |
| Zn 213.857              | 0.004                        |
| Y 371.029 (I/S)         | -                            |

## Standard and Sample Preparation

The calibration standards, blanks, and samples were prepared in accordance with the protocols outlined in ASTM Method D7111. Working standards were prepared by diluting oil-based stock standards (VHG, LGC Standards, Manchester, NH, USA) in PremiSolv™ (Conostan, AnalytiChem, Bale-D'Urfe, QC, Canada) prepared by weight. The Yttrium (Y) internal standard was added at a final concentration of 3 mg/Kg (ppm). A blank and a 2 mg/Kg mixed standard were prepared for all the listed analytes. Additionally, a 5 mg/Kg Phosphorus (P) standard was also prepared to provide a higher P calibration range.

Three middle distillate samples were sourced locally. These included a diesel fuel, an off-road (agricultural) diesel fuel, and Jet A fuel sample. Each of the samples were also prepared with a known spike (1 mg/kg, except P spiked at 2 mg/Kg) to evaluate system performance. The Y internal standard was also added to each of these fuel samples.

## Results

The MICAP was calibrated as described above and all wavelengths were viewed in the Radom Instruments Software (RIS) ProfilesView to verify no spectral interferences existed and to set peak integration points. The emission lines utilized for quantitation were selected from numerous line options available. The simultaneous, full wavelength range data collection also provided the ability during method development to confirm accurate measurements across multiple emission lines for each analyte element.

Figure 2 provides example emission peaks for the analytes of Mg, Cu, and Ti. Note the clear baseline resolution between the Ti 334.903 and 334.940 nm emission lines, demonstrating the high resolving power of the spectrometer. This capability helps ensure the analyte peaks of interest are not impacted by nearby spectral interferences from other potential trace elemental components.

The middle distillate fuel results and analytical spike results are displayed in Table 3. The three samples show very little elemental contamination is present. The sample spike results show excellent recoveries, with the exception of Boron (B) - for which stability / volatility issues are often quite common.

Internal standard recovery values were consistently within  $\pm 15\%$ , indicating that the system response was very stable between the calibration standards and the samples.




Figure 2. Profiles view of Mg, Cu & Ti emission lines in 1 mg/Kg spiked Jet A fuel. Sample emission lines displayed in brown with blank baseline indicated in green.

Table 3. Middle Distillate Fuel Analysis - sample and spike results (1 mg/Kg, P at 2 mg/Kg)

| Analyte Wavelength (nm) | Diesel Fuel   |                    | Agricultural Diesel Fuel |                    | Jet A Fuel    |                    |
|-------------------------|---------------|--------------------|--------------------------|--------------------|---------------|--------------------|
|                         | Conc. (mg/Kg) | Spike Recovery (%) | Conc. (mg/Kg)            | Spike Recovery (%) | Conc. (mg/Kg) | Spike Recovery (%) |
| Ag 328.068              | 0.038         | 100                | 0.010                    | 104                | -0.004        | 97                 |
| Al 396.152              | 0.028         | 105                | 0.028                    | 108                | -0.041        | 102                |
| B 249.772*              | 0.279         | 72                 | 0.067                    | 101                | 0.062         | 84                 |
| Ba 585.368              | 0.119         | 101                | 0.078                    | 105                | -0.308        | 102                |
| Ca 393.366              | 0.006         | 111                | 0.008                    | 114                | 0.015         | 110                |
| Cd 228.802              | 0.012         | 100                | 0.008                    | 102                | -0.021        | 97                 |
| Cr 428.973              | 0.003         | 102                | 0.005                    | 104                | 0.020         | 99                 |
| Cu 324.754              | -0.001        | 101                | -0.001                   | 104                | 0.000         | 98                 |
| Fe 259.940              | 0.007         | 102                | 0.042                    | 105                | -0.004        | 104                |
| K 766.490               | -0.077        | 88                 | -0.044                   | 87                 | 0.116         | 95                 |
| Mg 279.553              | 0.001         | 102                | 0.003                    | 104                | 0.002         | 103                |
| Mn 257.610              | 0.011         | 101                | 0.010                    | 103                | -0.010        | 103                |
| Mo 317.034              | -0.014        | 102                | -0.011                   | 104                | 0.015         | 97                 |
| Na 588.995              | 0.018         | 95                 | 0.012                    | 97                 | -0.015        | 95                 |
| Ni 300.249              | -0.004        | 106                | -0.007                   | 108                | 0.001         | 103                |
| P 213.618               | -0.076        | 106                | 1.646                    | 115                | 0.123         | 88                 |
| Pb 283.305              | 0.031         | 106                | 0.056                    | 104                | 0.017         | 96                 |
| Si 251.611              | -0.002        | 100                | 0.062                    | 99                 | 0.024         | 97                 |
| Sn 283.998              | 0.021         | 104                | 0.027                    | 106                | 0.008         | 100                |
| Ti 334.940              | -0.010        | 102                | -0.010                   | 104                | 0.001         | 104                |
| V 309.310               | -0.010        | 103                | -0.007                   | 105                | 0.017         | 103                |
| Zn 213.857              | 0.033         | 101                | 0.030                    | 103                | 0.012         | 100                |

\* Boron demonstrated carryover at low levels, impacting spike recoveries

## Conclusions

The MICAP-OES 1000 N<sub>2</sub> ICP-OES demonstrated its capability for the direct analysis of middle distillate fuels in accordance with ASTM Standard Method D7111. No carbon deposition from the samples was observed on the MICAP torch throughout this application work, demonstrating a significant impact on productivity as no user maintenance is required to clean the torch and then restart the system. The excellent accuracy and this stability confirms MICAP's capability for routine analysis of these fuel samples.

## References

1. ASTM D7111, Standard Test Method for Determination of Trace Elements in Middle Distillate Fuels by Inductively Coupled Plasma Atomic Emission Spectrometry (ICP-AES), ASTM Intn, West Conshohocken, PA, 2025, [www.astm.org](http://www.astm.org)
2. ASTM D5708, Standard Test Method for Determination of Nickel, Vanadium, and Iron in Crude Oils and Residual Fuels by Inductively Coupled Plasma (ICP) Atomic Emission Spectrometry, ASTM Intn, West Conshohocken, PA, 2025, [www.astm.org](http://www.astm.org)

Scan the code to access our website for additional information and our resources

